

Table 1 Identification of Gram-negative bacilli using VITEK® 2 ID-GNB cards inoculated with culture fluids from positive blood culture bottles

	Isolate tested		Correctly Identified		VNGNB		Not Identified		Misidentified	
	No.	No.	%	No.	%	No.	%	No.	%	
<i>Enterobacteriaceae</i>										
<i>Escherichia coli</i>	38	38	100.0	0	0.0	0	0.0	0	0.0	0.0
<i>Enterobacter faecalis</i>	14	13	92.9	0	0.0	0	0.0	1	7.7	
<i>Enterobacter cloacae</i>	3	3	100.0	0	0.0	0	0.0	0	0.0	0.0
<i>Klebsiella pneumoniae</i>	46	44	95.7	0	0.0	0	0.0	2	4.5	
<i>Klebsiella oxytoca</i>	1	1	100.0	0	0.0	0	0.0	0	0.0	0.0
<i>Salmonella</i> spp.	9	9	100.0	0	0.0	0	0.0	0	0.0	0.0
<i>Aeromonas</i> spp.	2	2	100.0	0	0.0	0	0.0	0	0.0	0.0
Total	113	110	97.3	0	0.0	0	0.0	3	2.7	
<i>Non-fermenters</i>										
<i>Acinetobacter baumannii</i>	42	40	95.2	2	4.8	0	0.0	0	0.0	0.0
<i>Chryseobacterium indologenes</i>	1	1	100.0	0	0.0	0	0.0	0	0.0	0.0
<i>Pseudomonas aeruginosa</i>	10	8	80.0	2	20	0	0.0	0	0.0	0.0
Total	53	49	92.5	4	7.5	0	0.0	0	0.0	0.0
<i>All types</i>	166	159	95.8	4	2.4	0	0.0	3	1.8	

VNGNB – various non-fermenting Gram-negative bacilli.