

**Table 2 Identification of Gram-positive cocci using VITEK® 2 ID-GPC cards inoculated with culture fluids from positive blood culture bottles**

Bacteria	Isolates tested	Correctly Identified		Not Identified		Misidentified		
	No.	No.	%	No.	%	No.	%	
<i>Enterococci</i>								
<i>Enterococcus faecalis</i>	11	10	90.9	1	9.1	0	0.0	
<i>Enterococcus faecium</i>	3	3	100.0	0	0.0	0	0.0	
Total	14	13	92.9	1	7.1	0	0.0	
<i>Micrococcus spp.</i>	3	3	100.0	0	0.0	0	0.0	
<i>Staphylococci</i>								
<i>Staphylococcus aureus</i>	31	29	93.5	2	6.5	0	0.0	
<i>Staphylococcus epidermidis</i>	2	2	100.0	0	0.0	0	0.0	
<i>Staphylococcus haemolyticus</i>	4	4	100.0	0	0.0	0	0.0	
Coagulase-negative staphylococci	8	5	62.5	0	0.0	3	37.5	
Total	45	40	88.9	2	4.4	3	6.7	
<i>Streptococci</i>								
<i>Streptococcus agalactiae</i>	2	2	100.0	0	0.0	0	0.0	
Group D streptococci	4	3	75.0	1	25.0	0	0.0	
<i>Streptococcus pneumoniae</i>	6	5	83.3	1	16.7	0	0.0	
Total	12	10	83.3	2	16.7	0	0.0	
<i>All types</i>	74	66	89.2	5	6.8	3	4.1	